BiFlu-Jen LC is a light curing dental varnish containing an active micronized suspension of Calcium Fluoride and Bioactive Calcium Fluoride-Phosphate containing Glass. The combined effectiveness of two ingredients leads to both, an immediate as well as long-term effects. Its optimal viscosity allows it to be easily applied to form a thin films and so serve as an insulator of exposed root surfaces. The Calcium Fluoride and Bioactive Calcium Fluoride-Phosphate containing Glass in BiFlu-Jen LC allow for deep fluoridation and remineralization when applied to floor of cavities or as a form of prophylaxis against caries and a source of fluoride ions when applied topically to the teeth.

• INDICATIONS

- Caries prophylaxis;
- Fluoridation and remineralization against caries;
- Management of exposed root surfaces;
- Protection of root caries lesion;
- Local reduction of bacterial activity;
- Topical application of fluoride following routine scale and polish;
- Post operative management of acid etched enamel surface;
- Insulation of tooth surface following mechanical treatment;
- Cavity liner, protective physical barrier;
- Dentinal tubule sealant.

• DIRECTIONS FOR USE

Surfaces to be treated are cleaned and thoroughly dried with oil-free air. Do not overdry the dentin surface.

When applied to the tooth surfaces, simply drip one or two drops of BiFlu-Jen LC onto a mixing pad and smear a thin layer upon the tooth surface with micro brush or other available applicator. The varnish is allowed to be absorbed for twenty seconds. Avoid forming thick film of varnish as this will tend easily to peel away. Use oil-free airjet to disperse BiFlu-Jen LC over treated surface. Cure the applied layer with available dental light unit. Note, that high-energy LED photo-polymerizer requires less time for curing the same portion of material (10 sec. is enough, as a rule). Follow to the manufacturer's instruction.

Remove the oxygen-inhibited layer with cotton pellet, control the margins and remove the residuals if necessary.

With correct application and proper mouth hygiene the film remains in place for several days.

For lining of cavities drip one or two drops onto an application brush and then use it to smear a layer of varnish onto walls

of the prepared tooth cavity. The varnish is allowed to be absorbed for twenty seconds and then dispersed with oil-free air. Cure the applied layer with available dental light unit. As a result BiFlu-Jen LC forms the water impermeable film. The routine therapeutic procedure is followed to BiFlu-Jen LC placing.

• SUPPLEMENTARY INDICATIONS

The application of BiFlu-Jen LC can be repeated two or three times at intervals of 5 – 7 days. Surfaces should be dried and treated with BiFlu-Jen LC again.

• PRECAUTIONS

- BiFlu-Jen LC is not recommended in case of dental fluorosis.
- BiFlu-Jen LC should not be used if patient has a known allergy or sensitivity to methacrylates or fluorides.
- Do not use BiFlu-Jen LC immediately after teeth whitening as peroxide residuals within dentin and enamel can significantly reduce the strength of resin bond to tooth tissue.
- Eugenol-containing temporary cements may also inhibit the polymerization of BiFlu-Jen LC.

• WARRANTY

The manufacturer warrants the quality of manufactured products. The adverse events inflicted by violation of user manual, storage conditions and other events inflicted by non-stipulated usage of the material are not the subjects of warranty. The customer is responsible for determination of suitability of this product for user's application. Warranty conditions: the product does not comply with requirements declared by manufacturer. In this case the manufacturer replaces the defective material within warranty period.

• LIMITATION OF LIABILITY

The manufacturer's liability is limited by only cases stipulated by direct legislation of the country.

• STORAGE

- Store at temperature of 4 – 25 °C. Protect from direct light. Close tightly the cap after use.
- Shelf life – 2 years.
- DO NOT FREEZE!

• RECYCLING

Dispose of the medical device in accordance with local / regional / national / international legal requirements.

• PACKING

- 3 g syringe with delivery tips.

For professional use in dentistry only!

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BiFlu-Jen LC

Light-curing fluoride varnish &
dentin tube sealant

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