Endogil-TC

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Endogil-TC is an excellent lubricant that facilitates the reaming of root canals which provides high-quality expansion, cleaning, shaping and antiseptic processing of even hard-to reach and crooked root canals.

●● Properties of material

- It doesn't require manual preparation (a ready form);
- It is delivered to the root canals through endodontic tip in a usual way;
- It has chelating properties, promotes more careful mechanical treatment of the root canal manually or with the help of rotational tools. Deletes the smear layer;
- It has antimicrobial effect on microflora of the root canal.

●● Composition (100 grams of the gel contain)

- Tetracycline hydrochloride – 2.5g;
- Ascorbic acid – 5g;
- Bioinert Solvent and the Biopolymer Matrix, Surfactant, Modifiers of a rheology and Fragrances – up to 100 g.

●● Indications to application

- The gel is used for instrumental treatment of root canals together with irrigational solutions (sodium hypochloride, hydrogen peroxide).
- It is especially recommended: for infected root canals / in case of destructive processes in periapical tissues.

●● Method of application and doses

First visit:

Gel is delivered into the canal by using an endodontical plastic tip. Mechanical treatment of the root channel is carried out according to the following scheme:

- Irrigation of the root canal with 2,5 – 3% solution of sodium hypochloride alternates with introduction of gel to the root canal;
- Instrumental formation of the root canal is carried out after entering the gel. Endogil-TC due to its chelating properties cleans the smear layer and facilitates passing through the root canal;
- The dentist should repeat the described procedures until he reaches the required degree of the root canal processing
- Finishing irrigation by sodium hypochlorite solution;
- Very careful irrigation of the root canal by distilled water;
- Temporary obturation of the root canal with calcium hydroxide for 1 – 2 weeks

Second visit:

After mechanical removal of calcium hydroxide from the root canal should be done the following:

- Irrigation of root canal with 2,5 – 3% sodium hypochlorite solution which alternates with gel applying to the root canal;
- Irrigation of the root canal by distilled water;
- Another temporary obturation of the root canal by calcium hydroxide (in case of destructive processes) or permanent obturation of the root canal.

Attention!

In cases of highly infected root canals usage of Jen-Metrohecor gel in the root canal under a temporary cement filling during 1 – 3 days can be required.

●● Active components of the gel and their pharmacological action

Tetracycline hydrochloride

Tetracycline hydrochloride suppresses bacterial growth by inhibiting protein synthesis. It binds to the 30S ribosomal subunit, and thereby inhibit binding of transfer RNA to the messenger RNA-ribosome complex. As a result, addition of amino acids to the growing peptide chain is prevented. At the concentrations achieved clinically, this substance is bacteriostatic.

It is active in relation to gram-positive (*Staphylococcus* spp., including producing penicillinase; *Streptococcus* spp., including *Streptococcus pneumoniae*); Haemophilus influenzae, *Listeria* spp., *Bacillus anthracis*) and to gram-negative microorganisms (*Neisseria gonorrhoeae*, *Bordetella pertussis*, *Escherichia coli*, *Enterobacter* spp., *Klebsiella* spp., *Salmonella* spp., *Shigella* spp.), and also *Rickettsia* spp., *Chlamydia* spp., *Mycoplasma* spp., *Treponema* spp.

Are steady against effect of medicine: *Pseudomonas aeruginosa*, *Proteus* spp., *Serratia* spp., majority of strains of *Bacteroides* spp. and mushrooms, small viruses.

●● Contraindications

Hypersensitivity to the drug and antibiotics, fungus diseases, abnormal liver functions, diseases of kidneys, leukopenia, pregnancy, lactation (for the period of treatment it is necessary to stop breast feeding), children's age up to 8 years.

●● Collateral reactions

The product can cause irritation of eyes, skin and mucous membranes.

●● Overdose

The casual or deliberate swallowing of a large amount of the gel can cause the following side effects:

In the alimentary system: nausea, vomiting, anorexia, abdominal pain, diarrhea, constipation, dryness in a mouth, glossitis, tongue discoloration, staining of teeth, esophagitis, transitional rising of content of hepatic transaminases, alkaline phosphatase, bilirubin, creatinine in blood.

In CNS: fainting, headache.

In a hematopoiesis system: neutropenia, thrombocytopenia, hemolytic anemia.

Allergic reactions: dermal eruption, itch, eosinophilia, Quincke's edema.

Dermatologic reactions: photosensitization.

Treatment: gastric lavage, if necessary symptomatic therapy.

●● Interaction with other medicines

Medicines containing ions of metals (calcium hydroxide or copper salts), form inactive chelates with tetracycline and in this connection it is necessary to avoid their simultaneous usage.

In case of local application in recommended doses interaction of the gel with other medicines isn't revealed.

●● Warranty

The manufacturer warrants the quality of manufactured products. The adverse events inflicted by violation of user manual, storage conditions and other events inflicted by non-stipulated usage of the material are not the subjects of warranty. The customer is responsible for determination of suitability of this product for user's application. Warranty conditions: the product does not comply with requirements declared by manufacturer. In this case the manufacturer replaces the defective material within warranty period.

●● Limitation of liability

The manufacturer's liability is limited by only cases stipulated by direct legislation of the country.

●● Storage

- Inappropriate storage conditions will reduce terms of use and can lead to deterioration in properties of material. Keep the product away from direct sunlight and do not leave the material under light for a long time.
- Store the product in the refrigerator at temperatures of 4 – 8 °C (39.2 – 46.4°F). Do not freeze!
- A warranty period of storage of the material – 12 months.

●● Recycling

Dispose of the medical device in accordance with local / regional / national / international legal requirements.

- Packing: The Endogil-TC – syringe 2 ml / Endodontic plastic cannulas (5 pieces).

For professional use in dentistry only!