form plast

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Zhermack
Dental

form plast

self-curing acrylic resin for patterns

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CE 0120

form plast

Self-curing acrylic resin for patterns

<p>Product characteristics: Self-curing acrylic resin for patterns. Burns without residues.</p>	
<p>Designation</p> <p>Supporting works in dental laboratories:</p> <ul style="list-style-type: none">- modeling of secondary parts of telescopic crowns, conical crowns and others, - modeling of inserts (inlay, onlay), tongue arch, palatine arch, clasps and other elements of skeletal prosthesis, - modeling of Maryland type adhesive bridges, - modeling of crown-root inserts, - structure modeling in post-implantation prosthesis, - temporary connections in soldering.	
<p>The works performed in the dental clinics:</p> <ul style="list-style-type: none">- modeling of endodontic posts and cores directly in the patient’s mouth	
<p>Contraindications: In people sensitive to the ingredients of the product allergic reaction may occur. If you notice any symptoms of allergic should stop using the product and consult a doctor.</p>	
<p>Assortment</p> <p>Set:</p> <p> [REF] V220Z01 form plast powder 30g + liquid 2x12ml</p> <p>Liquid:</p> <p> [REF] V220L02 form plast liquid 12ml</p>	
<p>Chemical composition</p> <p>Powder: acrylic polymers, benzoyl peroxide, pigments</p> <p>Liquid: methyl methacrylate, ethyl methacrylate, ethylene glycol dimethacrylate, N,N-dimethyl-p-toluidyne</p>	
<p>Instructions for use</p> <p>Procedure in technical and dental laboratories:</p> <p>Preparation: Class III and IV gypsum shall be used for working models, Stodent III and Stodent IV are particularly recommended. The gypsum model shall be isolated with compound which does not leave film on surface (DO NOT USE ISOLATORS BASED ON ALGINATES). Prepare two containers (glass, polyethylene (PE), ceramic containers or silicone cups). Transfer small amounts of powder and liquid to the prepared containers.</p> <p>Procedure: Material shall be applied on model with pencil brush. Soak the tip of pencil brush in liquid and then in small amount of powder. Place acrylic drop prepared in such manner on model. The quick gelation process causes that material does not flow from the place where it has been located. Spread material on the modeled element with pencil brush. This procedure shall be repeated till the required shape and thickness of the formed element is achieved. Material is hardened automatically in ambient temperature ca 23°C after ca 3-5 min. In case of modeling on metal (e.g. telescopic crowns), material shall be applied directly on the polished part, without any isolators. The pencil brush shall be washed in liquid, when necessary. The pencil brush shall be very carefully wash with liquid (methyl matacrylic) after completion of work.</p>	
<p>ATTENTION: Do not soak pencil brush with powder in bottle with liquid. The introduction of powder grains into liquid may cause not controlled polymerization process and the loss of properties by liquid.</p>	
<p>Processing and finishing: Carefully remove the formed mass from model (be careful not to damage the mould). Inspect visually the formed element and, if necessary, correct thickness and shape. Bring the formed element using the modeling wax to the final shape. Place the formed element on casting cone with wax, fasten in ring and then flood with protecting compound observing instructions of compound manufacturer. The formed element shall be cast and process in standard manner.</p>	
<p>Clinical procedure:</p> <p>Preparation: Prepare the root canals in the traditional way. Prepared root canal should remain wet, if necessary use wetting agents. Prepare a glass, polyethylene or other chemically resistant plate, HDPE pipette and brad plastic impressions burning without residues (Burn out post) to take "impression canal", two containers (glass, silicone or plastic) for a powder and liquid and brush modeling. Small amounts of powder and liquid move into prepared containers, and the prepared plate.</p>	
<p>Treatment:</p> <p>The brad plastic wetted in monomer dip in powder spilled on a plate and place in prepared root canal for about 10 seconds. Remove the post from the canal and moisten by monomer using a Pasteur pipette, and then dip in a powder, infiltrate the dry powder by monomer and place in the roote canal for approx. 10 seconds. Repeat the procedure for obtaining the correct impression of the root canal. Wait 4 minutes, cut off the sticking impression brad. Model the crown part post using a brush. Dip the tip of the brush in a liquid, and then in a small amount of powder. Put molded acrylic drop on modeled post. Fast gelling process causes that the material does not flow from the place where it was placed.</p>	
<p>Use a brush to spread the material on the modeled element. The procedure should be repeated until the desired shape and thickness of molded element is obtained.</p>	
<p>Brush should be cleaned and dried after each application. The hardening of the material goes on by itself at the temperature 23°C after approx. 3-5 minutes. Prepare the shaped element using cutters in order to obtain the proper shape and thickness. The finished restoration of crown-root pass to the dental laboratory.</p>	

Danger and counteraction

Information for dental technicians:

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: The place of product contact with skin wash thoroughly with soap with water.

P501 Dispose of contents/container to an approved waste disposal plant.

While working with the material do not: consume meals and smoke tobacco. Work should be carried out in good ventilated premises. Use stand ventilation.

Powder:

In case of contact of the powder particles with eyes the mechanical irritation may occur.

In case of contact acrylic powder with eyes wash it with water and contact physician.

Liquid:

Highly flammable liquid and vapor (containing methyl methacrylate), flash point +10°C. The place of liquid contact with skin wash thoroughly with soap water. In case of accidental swallowing contact the physician. Do not pour liquid into drainage system.

Storage and handling

Liquid is sensitive to influence of temperature and UV radiation under the influence of which the uncontrolled reaction of polymerization may occur.

Product components shall be stored in sales containers, in dark and good ventilated compartment, with temp. 5-25°C, place not accessible by children. Do not use the product after expiry date specified on the sales packaging and labels of the product ingredients

Disposal considerations

Dispose of in accordance with applicable regulations.

Carton packaging, labels, instructions - the possibility of recycling, dispose of as communal not dangerous waste.

Powder component and a packaging are not dangerous. Dispose of as communal not dangerous waste. In case of plastic container - the possibility of recycling.

Liquid and liquid packaging - hazardous waste. Pass to an authorized waste contractor.

NOTE!

Material for professional use only in technical and dental laboratories as well as dental clinics. Instructions in written form, oral or information relayed during Producers presentation result from actual level of dental technique and technology. It should be treated as obligatory information. Information’s obtained from third persons (due to their different level of authorisation) do not release the user from personal check whether the preparation is intended to the designated application. Due to the fact that the application takes place without the possibility of control from the producer’s side; therefore it limits the producers liability to the value of provided material.

We shall be very grateful for all comments related to products, please mail it to our address.

Zhermapol Sp. z o.o. has established Full Quality Assurance System certified to ISO 9001 / ISO 13485 and under surveillance by **Notified Body 0120, SGS United Kingdom Ltd. Systems and Services Certification.** Product has CE mark.

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