Jen-Line LCF is radiopaque light-curing one-component liner and base material containing calcium hydroxide, calcium hydroxyapatite, complex fluorides and adhesive promoters in a unique, urethane dimethacrylate formula. The liner has excellent bonding to dentine and is highly biocompatible.

•• INDICATIONS

- As a liner or base in prepared cavities.

•• DIRECTIONS FOR USE

Prepare the cavity, clean and dry it thoroughly. Care should be taken to ensure that no traces of eugenol are left in the cavity. Jen-Line LCF lining material should not be used in combination with the eugenol-containing materials. It is therefore recommended to use eugenol-free temporary cements.

Put the required portion of Jen-Line LCF on a mixing pad shortly before the application. Avoid exposure to daylight as it may cause premature setting of the material. Apply a layer of the Jen-Line LCF to the dentine using a small round condenser. Jen-Line LCF should be applied to dentine only. Do not apply it to enamel. The layer thickness should not exceed 1 mm. Should a thicker lining be necessary, apply and cure several subsequent layers. Each layer should be light-cured separately. Jen-Line LCF is cured by a special dental light unit (halogen, LED, etc.). The light guide should be positioned as close as possible to the liner. When the light guide is 2 mm above the liner surface, cure the liner for 20 seconds at least. When light guide is 5 mm above the liner surface, the curing time should be extended to 30 – 40 seconds. Due to high opacity of Jen-Line LCF the curing time should not be less than 20 seconds under any circumstances.

•• PRECAUTIONS

1. Jen-Line LCF has excellent biocompatibility. Nevertheless a possible allergic reaction cannot be excluded.
2. The liner should not be used for patients who are known to be allergic to derivatives of methacrylic acid.

•• SUPPLEMENTARY INDICATIONS

Jen-Line LCF is so strong after curing that a probe cannot be used to remove excess material as it would loosen the entire work. It is recommended to remove excess material from cavity linings and enamel using a rotary instrument.

•• WARRANTY

The manufacturer warrants the quality of manufactured products. The adverse events inflicted by violation of user manual, storage conditions and other events inflicted by non-stipulated usage of the material are not the subjects of warranty.

The customer is responsible for determination of suitability of this product for user's application. Warranty conditions: the product does not comply with requirements declared by manufacturer. In this case the manufacturer replaces the defective material within warranty period.

•• LIMITATION OF LIABILITY

The manufacturer's liability is limited by only cases stipulated by direct legislation of the country.

•• STORAGE

- Store at temperature of 4 – 25 °C. Protect from direct light. Close tightly the cap after use.
- Shelf life – 2 years.
- Cooling the material to 4 – 8 °C prolongs the shelf life.
- DO NOT FREEZE!

•• RECYCLING

Dispose of the medical device in accordance with local / regional / national / international legal requirements.

•• PACKING

- 3 g syringe with delivery tips.

JenD[®]

For professional use in dentistry only!

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