

pluser

LASER, FASTER, BETTER.



pluser




DON'T JUST STAY UP TO DATE. SET THE PACE.

Easier, faster, simply better. It's **Pluser**: the erbium laser that every modern dental practice should offer its clients, to **change their idea of a visit to the dentist's forever**. As time goes by, it is important to keep up with the clinical instruments designed to make your practice take speed, by maximizing benefits and minimizing worries. **Pluser is the answer that anticipates the question.**

LASER, FASTER, BETTER.

Pluser is a high-tech medical device that will launch the traditional dental practice into the future. The word Pluser means offering many things to your patients: **less pain, faster treatment, professional skill, smart versatility and continuous innovation in the treatment of soft and hard tissue, bone surgery and dermatology.**



EXTRAORDINARY IN YOUR DAILY ROUTINE.

High quality and top performance, from the most advanced surgical procedure to everyday treatments. **Pluser will shed new light** on commonly performed routine operations such as **tooth decay**, with a more relaxed and confident approach.

The erbium laser's high power speeds up cavity preparation and eliminates annoying vibration. The right amount of laser power will avoid any fractures that could weaken the structure of the tooth and preserve all healthy tissue.



MANY TREATMENTS, ONE **PLUSER.**



HARD TISSUE

- Caries removal
- Class I-V cavity preparation
- Etching
- Composite removal
- Sealing grooves
- Desensitization



SOFT TISSUE

- Surgery
- Frenectomy
- Implant uncovering
- Peri-implantitis
- Contouring
- Depigmentation



BONE SURGERY

- Split-crest
- Sinus lift
- Implant site preparation
- Crown lengthening
- Apicoectomy
- Bone implant



DERMATOLOGY

- Dermal surgery
- Non-ablative rejuvenation
- Fractional resurfacing

ANDROID ON BOARD.

Pluser uses the new proXimo software created by Doctor Smile on Android. You will be able to scroll and select the preset treatments with a simple swipe, thanks to the practical sliding menus. The 10.1" touch screen display lets you view every important operation parameter at all times. To get the most out of your Pluser, you can change the preset values with a simple touch.

- proXimo software
- Sliding menus
- HD 10.1" touch screen
- Four position display





ALWAYS CONNECTED.

Thanks to the Android platform and the **Wi-Fi connection**, Pluser is always on line, to offer you maximum support. The clinical protocols will be constantly updated with the latest research and innovation discovered by our laser tutors. The multimedia encyclopaedia offers a wide array of clinical videos. In order to solve any problem rapidly, Doctor Smile's service will be able to run a system diagnostic check remotely.

- **Direct Wi-Fi**
- **Live clinical protocols**
- **Tutorial videos**
- **24h service assistance**





RESEARCH AT ITS BEST.

Doctor Smile's Research and Development department has created new solutions to improve the way you work with Pluser. The innovative **MagicARM** supports the optic fibre so that you can operate in the most comfortable position, without any strain on your back. When not in use, the **MagicARM** fibre can be stored without taking up space. Pluser's original design will fascinate your patients, not scare them! Pluser can reach 12W of power, with a variable pulse frequency between 5 and 100 Hz, thanks to a new generation of erbium laser source. It is an extraordinary device, which performs exceptionally in any clinical situation, without limits.


- **MagicARM solution**
- **Flexible optic fibre delivery system**
- **12W peak power**
- **5 to 100 Hz pulse frequency**

THE SECRET IS IN THE PULSE.

Pluser's laser source was designed to create two different pulses: **Gaussian** and **Square**. The Gaussian pulse is ideal for dental applications because it guarantees tissue preservation and is well tolerated by patients: the energy is distributed gradually in order to reduce stress and pain.

The Square pulse is especially useful in dermatology because it can progressively release heat to coagulate and stimulate the dermis.

- **Gaussian Pulse for dental applications**
- **Square Pulse for dermal applications**
- **Complete Pulse shape personalization**
- **DEEP pulse for fractional resurfacing**



MULTITASKING IN A CLICK.

Pluser is a multi-disciplinary platform with numerous applications in dentistry and dermatology. The wide range of handpieces will accompany you in every treatment, with simplicity and effectiveness.

The **easy-click system** lets you switch handpieces easily, without wasting time. A quick touch lets you select your discipline: dentistry or dermatology. The graphic interface suggests which handpiece to use within the treatment list.

- **Easy-click system**
- **Dental handpieces: BOOST, Straight and 90°**
- **Dermal handpieces: full field and fractional**
- **Fractional scanner for skin resurfacing**



THE 10 TEMPTATIONS OF PLUSER.

- 1** **Fast and precise cutting:** The BOOST freebeam handpiece technology does not require tip changes: it operates without contact therefore reducing the risk of cross contamination.
 - 2** **Double impulse:** Pluser is the only laser that lets you choose the pulse shape: the Gaussian pulse for comfort in dentistry and the square pulse for effective treatments in dermatology.
 - 3** **Complete set of tools:** Pluser is a versatile, modular platform that allows fast handpiece and tip changes for different advanced applications, thanks to the easy-click connection system.
 - 4** **Built-in video guide:** The colour touch screen gives fast access to numerous preset treatments and clinical videos of many treatment procedures to help guide you step by step.
 - 5** **From dentistry to dermatology:** Pluser is a cutting edge tool for dermal surgery and skin rejuvenation thanks to its specific handpieces and treatment protocols.
 - 6** **Scanner-ready:** Pluser is designed to be fitted directly with a fractional scanner. It is the ideal tool for skin resurfacing, with excellent results and faster healing.
 - 7** **Optic fibre:** Flexible, practical and extremely durable: Pluser's optic fibre delivery system lets you operate freely, with greater comfort and efficiency.
 - 8** **Wireless footswitch:** Avoid messy cables in your work place: the wireless pedal allows you to power laser emission without wires.
 - 9** **Wi-Fi direct:** Pluser's Wi-Fi connection allows software upgrades, content upload and technical service, at any time.
 - 10** **100% Made in Italy:** Pluser is a CE certified medical device, designed and manufactured entirely in Italy.
-

DISCOVER THE **FUTURE** IN A SIMPLE **TOUCH.**

Pluser's smart: its intuitive, colour touchscreen gives fast and easy access to the preset treatments under the categories gum, tooth and bone.

Pluser is a **medical encyclopaedia**, always ready for consultation.

Its database arises from the lengthy experience of doctors and researchers on the use of the erbium laser.



Pluser has:

- A tablet interface
- Adjustable laser pulse
- Individualized parameter settings
- Wide range of handpieces
- Certified protocols
- Video and text content
- Animated simulations
- Clear clinical procedures
- Software upgrade for dermatology



PLUSER DERMA QUALITY IS NEVER ONLY SKIN DEEP.

Pluser gives you the opportunity to broaden your professional horizons, with the singular advantage of a complete treatment list. Pluser is so precise and versatile that it can be used in dermatology. With its specific tools your dental practice becomes a **multi speciality medical clinic that also offers aesthetic medicine, dermal surgery and fractional resurfacing.**

The Derma version includes a full field handpiece for skin lesions and non ablative rejuvenation.

The Fractional model contains the fractional handpiece and fractional scanner for skin resurfacing. What is the meaning of all this?

New treatments, new clinical cases, new patients.

pluser

LASER, FASTER, BETTER.





Smoothly shaped, easy to handle, limited size. This is Pluser. The **best laser technology** cased inside real Italian design. Front and rear handles, large wheels, limited weight: Pluser will be by your side during your greatest moments.

**DESIGN,
STYLE,
ATTITUDE.**

GREATER KNOWLEDGE MEANS EVOLUTION.



Do you want to start using laser technology but lack the experience? Would you like to learn more about a specific laser application? Doctor Smile, in collaboration with Master Class Academy, offers courses at various levels held by our skilled laser tutors that include the analysis of different clinical cases and hands on practical sessions. More information and the course agenda are available at masterclassacademy.it



Made in Italy



Pluser – Laser, faster, better. is a CE medical device entirely made in Italy.

Model	Doctor Smile Pluser
Code	L A ERT 001.1
Wavelength	2940 nm
Max power	12W
Max energy	600 mJ
Frequency	from 5 to 100 Hz
Size (WxHxD)	277x895x547 mm
Weight	35 kg
Medical /Laser class	II B / 4



All LAMBDA S.p.A. medical laser are certified with CE trademark in accordance with EEC Directives for medical device: EEC Directives 93/42

A M DP 184.1

Dealer



MORE INFO

pluserlaser.com • doctor-smile.com • info@doctor-smile.com

**DOCTOR SMILE IS
A TRADEMARK OF
LAMBDA SPA**

Via dell'Impresa, 1
36040 BRENDOLO (VI) Italy
T +39 0444.349165
F +39 0444.349954
info@lambdaspa.com
lambdaspa.com