

KENDA COMPLETE

2-STEP POLISHING SYSTEM

- Multifunctional system balanced for a wide range of restorative materials
- Economic polishers with high performance abrasives
- Available in the most common shapes
- Higher gloss performance compared to other non-diamond systems



SET REF 0500.012A



REF	reorder	0505	0508	0509	0510	0555	0558	0559	0560
ISO	1 - 658 204 ... 2 - 658 204 ...	030 523 060	243 523 050	303 523 100	030 523 090	030 513 060	243 513 050	303 513 100	030 513 090
	∅ × L mm	6,5 x 10	5 x 10	10 x 1	9 x 8,5	6,5 x 10	5 x 10	10 x 1	9 x 8,5
	μ grit	medium	light ruby colored	7500 - 10 000 RPM	10 000 RPM	fine	steel-grey	5000 - 7500 RPM	10 000 RPM
	color								
	rec.								
	max.								
	Refill Box of:	N Pressure:		Autoclavable:		Shank:		Medical device:	
	6, 12, 25, 100	0,3 - 2N		135°C, 275°F		RA, stainless steel		Class IIa	

